



**COSMOSIL**

HPLC Column for Separation of Soluble Carbon Nanotubes

# COSMOSIL CNT

- **Size-based separation of soluble carbon nanotubes**
- **Hydrophilic group coating silica packing material**
- **Three types of pore size (300 Å , 1000 Å , 2000 Å )**
- **High durability**





Image courtesy of Dr. Y. Ito, Univ. of Oxford

	CNT-300	CNT-1000	CNT-2000
Silica Gel	High Purity Porous Spherical Silica		
Average Particle Size	5 µm		
Average Pore Size	approx. 300 Å	approx. 1000 Å	approx. 2000 Å
Stationary Phase	Hydrophilic Group (Neutral)		
pH range	2 - 7.5		
Pressure	15 MPa and below		

COSMOSIL CNT series are ideal for separation of soluble carbon nanotubes based on sizes. COSMOSIL CNT series are packed with hydrophilic group-bonded silica packing material. The columns are specially designed to avoid adsorption of carbon nanotubes to silica support and thus ensure high resolution and maximum recovery of carbon nanotubes. COSMOSIL CNT series are available in three different pore sizes, 300 Å , 1000 Å and 2000 Å , respectively.

## Size-based Separation of Carbon Nanotubes

COSMOSIL CNT columns offered improved separation for DNA wrapped carbon nanotubes by connecting three columns with different pore sizes.




Column           COSMOSIL CNT-300 7.5 mm I.D. x 300 mm  
                       COSMOSIL CNT-1000 7.5 mm I.D. x 300 mm  
                       COSMOSIL CNT-2000 7.5 mm I.D. x 300 mm  
 Mobile phase    0.5 mmol/l EDTA, 200 mmol/l NaCl, 40 mmol/l Tris-Tris HCl(pH7.0)  
 Flow rate        1.0 ml/min  
 Temperature    40 °C  
 Detection        UV 260 nm(x 1), 350 nm(x 20)  
 Sample          DNA-wrapped CNT(CNT-DNA)  
 Injection vol.   500 µl  
 Pressure         8.1 MPa

Sample courtesy of Prof. H. Shinohara, Department of Chemistry & Institute for Advanced Research, Nagoya University

## AFM (Atomic Force Microscope) Analysis

The data below show analysis by AFM of CNT length of each fraction from the sample on the previous page.



## Ordering Information

Product Name	Column Size	Product Number
COSMOSIL CNT-300 Guard Column	7.5 mm I.D. x 50 mm	09194-81
COSMOSIL CNT-300 Packed Column	7.5 mm I.D. x 300 mm	09195-71
COSMOSIL CNT-1000 Guard Column	7.5 mm I.D. x 50 mm	09196-61
COSMOSIL CNT-1000 Packed Column	7.5 mm I.D. x 300 mm	09197-51
COSMOSIL CNT-2000 Guard Column	7.5 mm I.D. x 50 mm	09198-41
COSMOSIL CNT-2000 Packed Column	7.5 mm I.D. x 300 mm	09199-31

Other size may be available. Please enquire.

For research use only, not intended for diagnostic or drug use.