



Comprehensive two-dimensional GCxGC with conventional Inner Diameter columns: new column coated with two in series different stationary phases in a single fused silica tubing



Chiara Cordero¹; Carlo Bicchi¹; Patrizia Rubiolo¹; Mario Galli²; Stefano Galli²

¹ Dipartimento di Scienza e Tecnologia del Farmaco, Università degli Studi di Torino, Via P. Giuria 9, 10125 Torino, Italy
² MEGA, Via Plinio 29, I-20025 Legnano (Milano) Italy

Aims and scope

The optimization of a GCxGC system requires a complex approach since the separation in both dimensions is differently and independently influenced by column dimensions and stationary phases, temperature programming and carrier gas flow rates. In GCxGC, the columns of a set are combined in series and the possibility to apply a dose-to-optimal flow regime in both GCxGC dimensions by columns with a conventional ID was recently discussed [1] and the system performance evaluated through conventional chromatographic parameters (i.e. Peak Capacity (n) and Separation Measure (S_1 , S_2 and S_{GCxGC})) and a parameter specific to GCxGC such as orthogonality. The combination of homologous diameter columns, differing in stationary phase and film thickness, was shown to enable each chromatographic dimension to work with a dose to the optimal flow regime resulting in an improved phase selectivity (confirmed by system orthogonality estimation) that partially compensated the loss of efficiency due to wider 2D ID.

A new open tubular capillary column called **DN-UNIQUE™** or **MEGA-2D™** [2] coated in series with different film thickness of two different stationary phases in a single fused silica tubing is here described. This column is here shown to improve GCxGC performance since it avoids the use of unions (press-fits or low dead volume connections) between the first and the second dimension thus eliminating a possible source of leaks and reducing band broadening effects.

The advantages of the single column are here shown through the results of a set of samples including n -alkanes (C9-C25), Fatty Acids Methyl Esters (FAME, C4:0-C24:0), hydrocarbons (Boiling Point range 36-254°C), and suspected volatile allergens.

References:
[1] C. Cordero, C. Bicchi, E. Libertò, B. Sgorbini, S. Galli, M. Galli, P. Rubiolo, Proceedings of Dalan International Symposia and Exhibition on Chromatography, June 4-5, 2001, Dalan, China
[2] Patented System by DANI Instruments SpA.

Experimental

Samples

Pure standard samples of **n -alkanes** (from C9 to C25), and **suspected volatile allergens**: limonene, phenylacetaldehyde, linalool, benzyl alcohol, estragole, methyl 2-cyanoacrylate, citronellol, geraniol, citral, cinnamic aldehyde, hydroxycitronellal, anisyl alcohol, cinnamic alcohol, eugenol, methyl eugenol, α -isomethylionone, isoeugenol, butylphenyl methylpropional (lilial), coumarin, amyl cinnamic aldehyde, farnesol, amyl cinnamic alcohol, hydroxyisohexyl-3-cyclohexene carboxaldehyde (lyral), heptyl cinnamic aldehyde, benzyl benzoate, benzyl salicylate, benzyl cinnamate, and 1,4-dibromobenzene (1STD-1), 4,4'-dibromodiphenyl (1STD-2) were supplied by Sigma-Aldrich (Milan, Italy). Solvents (cyclohexane, n -hexane, acetone) were all HPLC-grade from Riedel-de Haen (Seelze, Germany). Standard stock solutions and standard working solutions were stored at -18°C .

The **Fatty Acids Methyl Esters mixture** was purchased from Supelco (Milan, Italy) and the hydrocarbons **"Quantitative Reference Standard 512"** mixture (Boiling Point range 36-254°C) was purchased from AC Analytical Controls (Rotterdam, The Netherlands) and consisted of: cyclopentane, n -pentane, cyclohexane, 2,2-dimethylbutane, n -hexane, 1-hexene, methylcyclohexane, 4-methyl-1-hexene, n -heptane, 1,2-dimethylcyclohexane, 2,2,4-trimethylpentane, n -octane, 1,2-trimethylcyclohexane, nonane, n -decane, n -undecane, n -dodecane, benzene, toluene, trans-decalin, n -tetradecane, ethylbenzene, o -xylene, n -propylbenzene, 1,2,4-trimethylbenzene, 1,2,3-trimethylbenzene, 1,2,4,5-tetramethylbenzene, pentamethylbenzene. Components range from a minimum of 0.1% to a maximum of 5% w/w. The complex **aroma** sample was supplied by Robertet S.A. Grasse (Cedex) France.

Instrumental set-up

Comprehensive GCxGC/qMS analyses were carried out on an Agilent 6890 GC coupled with a 5975 MS detector (Agilent, Little Falls, DE, USA) operating in E.I. mode at 70 eV. Ion source temperature: 230 $^\circ\text{C}$; Quadrupole temperature 150 $^\circ\text{C}$; Transfer line: 280 $^\circ\text{C}$. An automatic tuning was used. Scan range was from 35 m/z to 300 m/z with a scan rate of 10000 amu/s. The system was provided with a two-stage thermal modulator, **Figure 1** (KT 2004 loop modulator from Zoex Corporation, Lincoln, NE, USA) cooled with liquid nitrogen and with the hot jet pulse time set at 250 ms. Data acquisition was by Agilent - MSD Chem Station ver D.02.00.275 and data elaboration by GC-Image ver 1.8.866 LLC Lincoln (NE) USA.

GCxGC Operating conditions

Table 1 reports column characteristics and operative conditions adopted in this study. All **DN-UNIQUE™** or **MEGA-2D™**, **Figure 2**, columns were from MEGA (Legnano (Milan)-Italy).

One micro liter of each sample solution was automatically injected into the GC instrument by an Agilent ALS 7683B under the following conditions: injector: split/splitless in split mode, split ratio: 1/200, injector temperature: 280 $^\circ\text{C}$; Carrier gas: Helium.

Temperature programme: from 50 $^\circ\text{C}$ (1 min) to 280 $^\circ\text{C}$ (5 min) at 3 $^\circ\text{C}/\text{min}$. The modulation period was set at 4 s.

Results and Discussion

Peak Capacity and Separation Measure

Peak Capacity (n) and Separation Measure (S), were adopted to evaluate the separation power of the **DN-UNIQUE™** or **MEGA-2D™** using the average α values obtained with the separation of the n -alkanes test mixture.

Peak capacity (n), was defined by Giddings as the maximum number of peaks in a selected time interval separated with a given resolution. It was calculated through the following equation:

$$n = \Delta t/w_b \quad \text{Eq1}$$

where Δt is the time interval and w_b is the base peak width that can be assumed to be four times the standard deviation (σ) of the peak. The separation measure S introduced by Blumberg *et al* [1], is an additive quantity but it is representative of a separation time interval which is equal to the sum of the separation measures of non-overlapping α -wide subintervals that can be used with any shape of chromatographic peaks. It was calculated using the following equation:

$$S = \Delta t/\sigma_{av} \quad \text{Eq2}$$

where Δt is the arbitrary time interval limited by two peaks a and b , $\Delta t = t_b - t_a$, and σ_{av} is the average σ of a and b , $\sigma_{av} = (\sigma_a + \sigma_b)/2$. In 2003 Blumberg [2] extended the S concept to a GCxGC separation introducing the S_{GCxGC} that corresponds to the product of the separation measure of each chromatographic dimension.

The **Experimental S_{GCxGC}** was calculated on C9-C25 separation intervals respectively defined as: $^1D(t_{C25} - t_{C9})$ and $^2D(t_{C25} - t_{C9})$ using the $^1D \sigma_{av}$ was obtained from the analysis of n -alkanes with each column without modulation and with the other chromatographic conditions constant, $^2D \sigma_{av}$ was measured by analyzing n -alkanes with each column combination and a modulation period of 4 s.

Table 2 reports n values and **Experimental S_{GCxGC}** measured on the same time intervals and α parameters. **S_{GCxGC} values are perfectly comparable to those of "classical" GCxGC column setting (i.e. with a narrow bore column in the 2D) and compatible with a separation power suitable for complex samples.**

[1] L.M. Blumberg, M.S. Kline, J. Chromatogr. A, 933 (2002) 1
[2] L.M. Blumberg, J. Chromatogr. A, 985 (2003) 29

Table 2: 1D absolute retention time (t_b), $^1D \sigma$ ($^1D \sigma$), 2D absolute retention time (t_b), $^2D \sigma$ ($^2D \sigma$), 1D hold-up time (t_{H0}), 1D analysis time (t_{AN}), 1D Separation Measure (S_1), 2D Separation Measure (S_2), GCxGC Separation Measure (S_{GCxGC}), GCxGC Peak Capacity (n) estimated on the basis of C_9 and C_{25} elution time intervals for each set.

	Reference solute n -C ₉				Reference solute n -C ₂₅				Practical			
	t_b (s)	$^1D \sigma$ (s)	t_b (s)	$^1D \sigma$ (s)	t_b (s)	$^2D \sigma$ (s)	t_b (s)	$^2D \sigma$ (s)	S_1	S_2	S_{GCxGC}	n
Thick OV1701	398	0.94	0.59	0.04	3580	2.17	4.00	0.00	2043	84	171952	10749
Thick OV17	394	1.02	0.59	0.06	3598	2.25	4.00	0.00	1961	56	109556	6872
MEGA-2D™ OV1-OV1701	380	0.91	0.59	0.04	3579	2.09	4.00	0.00	2130	93	197086	12318
MEGA-2D™ OV1-OV17	383	1.00	0.59	0.04	3574	2.09	4.00	0.00	163	68	193979	8711

Separation Space Used and Peak Spreading

The % of usage of the separation space [1,2] was used to investigate the degree of correlation between the two dimensions on the basis of the peak distribution on the chromatographic plane. This parameter was proposed by Ryan *et al* [1] and is a practical measure of the degree of orthogonality, a fundamental aspect for a GCxGC separation. It measures the ratio between the area occupied by solute separation and the unused separation space beneath the 2D dead time. **Table 3** reports the amount of separation space used (referred to three different model test mixtures: suspected volatile allergens, hydrocarbon and FAME mixtures). Two series of data are reported: the first one, called **Separation space used***, includes the wrapped around peaks, i.e. all separated peaks are included, and the lower unused area of the retention plane is included within the least-retained peak. The second one excludes them. The net separation space through which data were normalized, was referred to 2D column dead time ($^2D t_0$) calculated by Poiseuille's law. **Experimental data clearly show that the % of usage of the separation space is maximized even if conventional i.d. columns are used. A fairly separation of the suspected allergens standard mixture is reported in Figure 3.**

[1] D. Ryan, P. Morrison, P. Marriott, J. Chromatogr. A 1071 (2005) 47
[2] C. Cordero, P. Rubiolo, B. Sgorbini, M. Galli, C. Bicchi, J. Chromatogr. A 1132 (2006) 26

	1D last eluted (min)	1D least eluted (s)	2D last eluted (s)	Total available separation space (s ²)	Separation space used*	Separation space used	% of usage
FAME test mixture							
Thick OV1701	44.67	0.59	0.29	3.95	1.07	0.99	99
Thick OV17	46.78	0.59	0.15	3.76	0.97	0.93	93
MEGA-2D™ OV1-OV1701	43.47	0.59	0.56	3.96	0.99	0.96	96
MEGA-2D™ OV1-OV17	45.00	0.59	0.48	3.88	0.97	0.96	96
Hydrocarbons test mixture							
Thick OV1701	22.45	0.59	0.24	3.95	4.593	1.09	0.99
Thick OV17	22.67	0.59	0.24	3.76	4.638	1.03	0.93
MEGA-2D™ OV1-OV1701	22.33	0.59	0.21	3.89	4.569	1.08	0.96
MEGA-2D™ OV1-OV17	22.93	0.59	0.08	3.96	4.692	1.14	0.99
Suspected allergens test mixture							
Thick OV1701	44.67	0.59	0.29	3.95	9139	1.07	0.99
Thick OV17	46.78	0.59	0.05	3.76	9571	1.09	0.93
MEGA-2D™ OV1-OV1701	43.47	0.59	0.56	3.96	8893	1.00	0.99
MEGA-2D™ OV1-OV17	45.00	0.59	0.48	3.88	9207	1.00	0.96

* included peaks beneath the 1D hold-up time

Table 3: Amount of separation space used and % of usage of the separation space available calculated for the allergens, the Fatty Acids Methyl Esters and hydrocarbons test mixtures.




Figure 2: Possible Points of Leaks and Handling Difficult





Figure 3: Contour plots of the Suspected Allergens standard mixture (3a), Hydrocarbons "Reference standard 512" (3b) and C:20 FAME cluster (3c) (methyl arachidate (C20:0), cis-11-eicosenoic acid methyl ester (C20:2w6), cis-6,11,14-eicosatrienoic acid methyl ester (C20:3w3), cis-11,14,17-eicosatrienoic acid methyl ester (C20:3w3), methyl arachidate (C20:4w6), cis-5,8,11,14,17-eicosapentaenoic acid methyl ester (C20:5w3) obtained with MEGA-2D™ OV1-OV1701

Conclusions

Experimental data demonstrate that the two **MEGA-2D™** or **DN-UNIQUE™** columns tested, with 0.25 mm homologous diameter coated with different film thickness of two stationary phases in series corresponding to the two chromatographic dimensions, operating with a flow regime close to the optimal linear velocity, gives high separation power and phase selectivity (confirmed by system orthogonality estimation). Moreover, a suitable tuning of the elution temperature in combination with the choice of a suitable thicker film in the 2D column makes it possible:

- to compensate the loss of separation efficiency, even because the increase in elution temperature results in a narrower 2D peak width and a higher peak capacity,
- to enhance the separation space used and
- to obtain a suitable number of modulated peaks enabling a reliable quantitation also for trace analytes. Last but not least, homologue diameter column combinations improve peak compatibility with quadrupole MS detection giving a wider 2D peak width, when compared to those obtained with a "conventional" narrow bore 2D column.

Acknowledgments

This research was carried out within the project entitled: "Sviluppo di metodologie innovative per l'analisi di prodotti agroalimentari" (FIRB Cod.: RBIP065XMR_002) of the Ministero dell'Istruzione, dell'Università e della Ricerca (MIUR) (Italy).



A new column coated with two in series different stationary phases in a single fused silica tubing of conventional inner diameter for Comprehensive Two-Dimensional GCxGC

Chiara Cordero¹; Carlo Bicchi¹; Barbara Sgorbini¹; Patrizia Rubiolo¹; Mario Galli²; Stefano Galli²
 1 Dipartimento di Scienza e Tecnologia del Farmaco, Università degli Studi di Torino, Via P. Giuria 9, 10125 Torino, Italy
 2 MEGA, Via Plinio 29, I-20025 Legnano (Milano) Italy

Aims and scope

Several parameters influence the optimization of a GCxGC system including column dimensions and stationary phases, temperature programming and carrier gas flow rates. The use of two conventional 1D columns in both GCxGC dimensions affording a close-to-optimal flow regime was recently discussed [1] by evaluating Peak Capacity (n), Separation Measure ($S1$, $S2$ and S_{GCxGC}) and orthogonality. This combination provided an improved phase selectivity compensating the loss of efficiency due to higher 2D ID.

A new open tubular capillary column [2], called **DN-UNIQUE™** or **MEGA-2D™** coated with two different stationary phases (differing in composition and film thickness) in a single fused silica tubing is here presented. This column is an effective GCxGC improvement since it avoids unions between the two dimensions thus eliminating possible leaks and reducing band broadening effects. The performance of this column is here compared to that of a conventional 1D 0.25 mm ID and 2D 0.10 mm ID column set up. **DN-UNIQUE™** or **MEGA-2D™** peak capacity, orthogonality, peak area reproducibility and linearity have been here evaluated through the analysis of two test mixtures (Fatty Acids Methyl Esters and volatile suspected allergens) and through a target analysis approach aimed at quantifying volatile suspected allergens in medium-complexity fragrances.

[1] C. Cordero, C. Bicchi, P. Rubiolo, M. Galli, S. Galli, 59th Pittsburgh Conference on Analytical Chemistry and Applied Spectroscopy, March 1-7, 2008 New Orleans, Louisiana, USA
 [2] Patent pending by DANI Instruments SpA and MEGA snc

Experimental

Samples

Pure standard samples of **volatile suspected allergens** were supplied by Sigma-Aldrich (Milan, Italy).
Solvents (cyclohexane, ethyl acetate) were all HPLC-grade from Riedel-de Haen (Seelze, Germany).
The Fatty Acids Methyl Esters mixture was purchased from Supelco (Milan, Italy).
The Fragrance sample was supplied by Robertet S.A. Grasse (Cedex) France.

Instrumental Set-up

Comprehensive GCxGC/MS analyses were carried out on an Agilent 6890 GC coupled with a 5975 MS detector (Agilent, Little Falls, DE, USA) operating in E.I. mode at 70 eV, ion source temperature: 230 °C, Quadrupole temperature 150 °C, Transfer line: 280 °C. An automatic tuning was used. Scan range was from 35 m/z to 300 m/z with a scan rate of 10000 amu/s. The system was provided with a two-stage thermal modulator (RT 2004 top modulator from Zoex Corporation, Lincoln, NE, USA) cooled with liquid nitrogen and with the hot jet pulse time set at 250 ms. Data acquisition was by Agilent - MSD Chem Station ver.02.00.275 and data elaboration by GC-Image ver. 1.8666 LLC Lincoln (NE) USA.



GCxGC Operating conditions

Table 1 reports column characteristics and operative conditions adopted in this study. All **DN-UNIQUE™** or **MEGA-2D™** Figure 1, columns were from MEGA (Legnano (Milan) Italy). One micro liter of each sample solution was automatically injected into the GC instrument by an Agilent ALS 7683B under the following conditions: injector: split/splitless in split mode, split ratio: 1/200 or 1/20, injector temperature: 280 °C; Carrier gas: Helium. Temperature programme: from 50 °C (1 min) to 280 °C (5 min) at 3 °C/min. The modulation period was set at 4 s.

Table 1

Column Set N°	First dimension column (length m x ID mm, df µm)	Second dimension column (length m x ID mm, df µm)	P ₀ (kPa)	u ₁ (cm s ⁻¹)	u ₂ (cm s ⁻¹)	t _D Hold-up time (s)
1	OV1-25x0.25, 0.15	OV1701-1.0x0.10, 0.10 (loop 1.0x0.10)	298	15	147	0.7
2	MEGA 2D™ OV1, 0.15	OV1701-1.0x0.10, 0.10 (loop 1.0x0.25)	175	48	149	1.7
3	MEGA 2D™ OV1, 0.15	OV1701-1.0x0.10, 0.10 (loop 1.0x0.25)	175	48	149	1.7

Figure 1



Results

Basic Performance: Peak Capacity

GCxGC net Peak Capacity (n_{GCxGC}) was adopted to evaluate the separation power of the **DN-UNIQUE™** or **MEGA-2D™**.

Peak capacity (n), defined by Giddings as the maximum number of peaks in a selected time interval separated with a given resolution, was calculated for each chromatographic dimension through the equation:

$$n = \Delta t / w_b$$

where Δt is the time interval, w_b is the base peak width assumed to be four times the standard deviation (σ) of the peak.

Net peak GCxGC capacity (n_{GCxGC}) is the product of n calculated for each chromatographic dimension. It was referred to the separation intervals between the first and the last eluted compound of the two test mixtures in the 1D, i.e. $2D_{GCxGC} = (t_{R, last} - t_{R, first}) / n_{GCxGC}$, and for the 2D it was comprised between the hold-up and the total analysis time, i.e. $2D_{GCxGC} = (t_{total} - t_{hold-up}) / n_{GCxGC}$. $2D_{GCxGC}$ was estimated by analyzing test mixtures without modulation and keeping the other chromatographic conditions constant while $2D_{GCxGC}$ with a modulation period of 4 s.


Table 2 reports 1D and 2D n calculated for the first and the last eluted components of the two test mixtures while Figure 2 reports n_{GCxGC} values calculated for the column set under study.

Table 2

Set N°	Benzyl alcohol		Benzyl Cinnamate		n_{GCxGC}
	1D (s)	2D (s)	1D (s)	2D (s)	
1	1.81	0.08	2.29	0.08	3214
2	1.53	0.20	2.32	0.18	952
3	1.53	0.16	2.32	0.17	1112

Set N°	Methyl hexanoate		Methyl tetraosanoate		n_{GCxGC}
	1D (s)	2D (s)	1D (s)	2D (s)	
1	2.04	0.03	2.30	0.03	11961
2	1.66	0.09	1.79	0.05	4581
3	1.66	0.08	1.79	0.04	5056

Figure 2



Basic Performance: Separation Space Used

The % of usage of the separation space [1,2] was used to investigate the degree of correlation between the two dimensions on the basis of the peak distribution on the chromatographic plane. This parameter is a practical measure of the degree of orthogonality and indicates the ratio between the area occupied by solute separation and the unused separation space beneath the 2D dead time. Figure 3 reports the amount of separation space used referred to volatile suspected allergens and FAME test mixtures. The net separation space through which data were normalized, was referred to 2D column hold-up time ($2D_{w_0}$).

Experimental data show that the % of usage of the separation space is maximized with conventional 1D columns, **DN-UNIQUE™** or **MEGA-2D™**, because of the improved exploitation of the 2D stationary phase. A fair separation of the suspected allergens standard mixture is reported in Figure 4.

[1] D. Ryan, P. Morrison, P. Marriott, J. Chromatogr. A 1071 (2005) 47
 [2] C. Cordero, P. Rubiolo, B. Sgorbini, M. Galli, C. Bicchi, J. Chromatogr. A 1132 (2006) 26

Figure 3





Figure 4



Target Analysis: Precision and Linearity

Volatile suspected allergen quantitative determination is part of the fragrance quality control assessment and should take into account the highly variable distribution of fragrance components and adopt experimental conditions to compensate peak distortions due to both column overloading and strong retention effects. Widely used orthogonal stationary phase combinations, such as 1D OV1/2D CW20M, were not effective enough in particular with the polar analytes affected by a strong 2D retention due to the inappropriate elution temperature of the 2D column. These effects are well overcome by the adoption of OV1701 or OV17 but the higher correlation (low orthogonality) between the two dimensions for the 0.25/0.10 mm ID column settings need for slower temperature rates, and consequently, higher analysis times. The exploitation of 2D stationary phase selectivity showed by 0.25 mm homologous 1D column combination, resulting in significant peak spreading over the chromatographic plane, and the higher 2D column loadability suggested to test **DN-UNIQUE™** or **MEGA-2D™** for an effective separation of target allergens and to evaluate if their peak capacity was high enough to separate target analytes in a medium complexity fragrance.


Correlation coefficients (R^2) estimated by regression analyses over a 50-2 mg/L range and precision results (expressed as RSD%) referred to the 2D area of each analysis measured over six replicates, are reported in Table 3.

Results show that **MEGA-2D™** OV1701 operating at controlled flow, temperature and modulation period conditions, can successfully be used for the target analysis of volatile suspected allergens in medium complexity fragrance (see Figure 5). **MEGA-2D™** is characterized by an improved exploitation of stationary phase selectivity and an increased 2D column loadability; its limited net peak capacity, if compared with an equivalent conventional column setting limits its use to samples with medium complexity.

Table 3

Component Name	SIM m/z ions	Regression analyses R ²		2D Area Precision RSD%	
		Set N° 1	MEGA-2D	Set N° 1	MEGA-2D
amylcinnamic aldehyde	202, 201, 129	0.992	0.996	2.05	2.03
anilyl alcohol	138, 137, 109	0.973	0.994	2.15	1.93
benzyl alcohol	108, 79, 107	1.000	0.997	2.61	0.52
benzyl benzoate	105, 212, 194	0.997	0.990	2.05	1.62
benzyl salicylate	91, 228, 65	0.999	0.994	1.64	0.93
cinnamic alcohol	92, 134, 115	0.980	0.999	1.37	2.60
cinnamic aldehyde	131, 132, 103	0.991	0.998	5.54	4.59
coumarin	146, 118, 89	0.996	0.998	2.05	3.21
farnesol isomer I	69, 93, 81	0.986	0.999	2.15	2.15
farnesol isomer II	69, 93, 81	0.988	0.999	2.61	2.06
geraniol	69, 123, 93	0.991	0.998	2.05	2.05
hexylcinnamic aldehyde	216, 215, 129	0.980	0.999	1.64	1.62
linalone	154, 149, 131	0.999	0.998	1.37	1.36
linalool	68, 93, 67	0.999	0.999	5.54	1.98
phenylacetalddehyde	91, 120, 92	0.995	0.997	1.77	3.52

Figure 5



Conclusions

Experimental data demonstrate that the two **MEGA-2D™** or **DN-UNIQUE™** columns tested, with 0.25 mm homologous diameter coated with different film thickness, with two stationary phases in series operating with a flow regime close to the optimal linear velocity, could successfully be used for specific applications. This unique column avoids unions between the two dimensions and, as a consequence, possible leaks and band broadening effects. Their performances were verified through the definition of several method performance parameters in a target analysis approach for several allergens quantification in medium complexity fragrances. Linearity over the working range and precision were, above all, issues of interest even if compared with GCxGC separations performed with conventional column combinations. Despite their reduced peak capacity the proper exploitation of 2D stationary phase can compensate this loss giving reliable results in shorter analysis times.

Acknowledgments

This research was carried out within the project entitled: "Sviluppo di metodologie innovative per l'analisi di prodotti agroalimentari" (FIRB Cod.: RBIP06SXMR_002) of the Ministero dell'Istruzione, dell'Università e della Ricerca (MIUR) (Italy).