

CTS-10988

2690/2695 Performance Maintenance Kit
Comparable to: Waters® WAT270944



Technical Data Sheet

1.0 Description

This PM kit contains all the parts necessary to keep Waters 2690, 2690D, 2695, 2695D and Alliance HPLC systems working at peak performance. It is comparable to Waters part number WAT270944.

2.0 Specifications

2.1 Components

Qty	Part Description
2	Sapphire Plunger (CTS-10476)
1	Seal Wash Plunger Seal Kit, 2/pk (CTS-10560)
1	Head Plunger Seal Kit, 2/pk (CTS-10566)
1	Wash Tube Seal Kit, 4/pk (CTS-10597)
1	Sparge Diffuser (CTS-10984)
1	Face Seal, 4/pk (CTS-10989)
4	Solvent Reservoir Filter, 10µm (CTS-5531)
1	250µL Syringe (CTS-3109)
1	2µm Precolumn Filter Insert (CTS-10919)
2	Seal Wash Face Seal (CTS-10532)
1	PP Check Valve Cartridge, 2/pk (CTS-10600)
1	Seal Pack Rebuild Kit (CTS-10839)

2.2 Wetted Materials

18-8 SS
302 SS
316 SS
316 SST



Borosilicate Glass
Delrin®
Elgiloy Spring
HDPE
MP35N SS
PCTFE
PEEK™
PTFE
Ruby
Sapphire

3.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.

4.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at sales@sciencix.com to report the issue. We will assist you with a resolution in a timely fashion.