



Q-Sonic Ultrasonic Cleaner Solution

QLA's NEW Ultrasonic Cleaning Solution was designed specifically for your busy lab environment. **There is no longer a need to change cleaning solutions for different materials!** Our proprietary, concentrated formula effectively cleans 316 Stainless Steel, PEEK, PTFE, silicone, and other common lab accessory materials. Q-Sonic Ultrasonic Cleaner Solution is most effective when used in our frequency and temperature optimized Ultrasonic Cleaner but can be used with any sonicating water bath.



- ✓ **All-in-one formula saves you money.**
- ✓ **Proven more effective than DI Water alone.**
- ✓ **Non-Hazardous.**
- ✓ **No special disposal necessary.**
- ✓ **Non-Flammable.**
- ✓ **No staining or discoloration.**
- ✓ **No odor.**
- ✓ **Non-Oxidizing.**

QLA Part# QSONIC Q-Sonic Ultrasonic Cleaner Solution, 8oz Bottle

QLA Part# QSONIC-3 Q-Sonic Ultrasonic Cleaner Solution, 8oz Bottle (3-Pack)

YOUR DISSOLUTION SPECIALISTS